


PREMIUM WOVEN SURGICAL ABDOMINAL BINDERS

Sizes: Small - Large
Color: White

Three-Panel, 9" Model No: 34-630



Four-Panel, 12" Model No: 34-640



The Premium Woven Abdominal Binders provide support and compression to abdominal muscles weakened by strain without restricting breathing or muscle redevelopment. Abdominal binders enhance compliance following open abdominal surgical procedures and can help prevent sedentary complications such as pneumonia. Wearing a binder encourages ambulation earlier and allows for more comfortable and controlled breathing. Higher 12" 4-panel binder for taller patients or long-line surgeries. Easy and inconspicuous to wear under most clothing. For Men and Women.

FEATURES:

- **True Woven Elastic** – Gives firm and consistent compression power and can be cut without unraveling to accommodate drainage tubes.
- **Paneled Hinged Stitch** – Creates separate panels for a more conforming, tapered fit with flexibility. Does not pinch the skin.
- **Hook Compatible Exterior** – Allows for a wider range of adjustability and easy application. Holds up through repeated use.
- **Latex Free and Plush Lining** – Comfortable and non-irritating for patients and caregivers with sensitivity to latex rubber.

HOW TO APPLY:

Individual Application: Apply in a standing or laying position. Open belt to its fullest length with the care tag on the inside. Center the height of the abdominal binder at the navel with the seam centered in the back. The hook should be on the right side. Stretch both sides and wrap the sides without the hook around the body first. Secure the hook onto the support where needed. Adjust as necessary.

Caregiver Application: With patient lying flat on his/her back in bed, open the binder to its fullest length, positioning it underneath the patient with the plush side against the patient's skin. Center the height of the abdominal binder at the navel with the seam centered in the back. The hook should be on the right side of the patient. Stretch the left side of the binder and pull it across the patient's ribs/abdomen, using the weight of the patient to hold it in place. Holding the left side in place on the patient's ribs/abdomen, stretch the right side the same way securing to the loop. Adjust as necessary.

CARE INSTRUCTIONS:

Fasten closure and hand wash in cold water with mild soap. Wrap in towel to remove excess water (do not twist or wring). Air dry away from heat. Do not tumble dry. Do not bleach.

SIZING:

Measure around the fullest part of the abdomen.

Size	Fits
Small	30 - 45"
Medium	46 - 62"
Large	60 - 75"



800-327-4110

Website: www.flaorthopedics.com

©2004 FLA ORTHOPEDICS, INC.